## How to use the NOVASINA accessory “Sorption Isotherm Set”

### **Water in Foods**

Most materials and all food products contain water: bound water and free water.

The WATER or MOISTURE CONTENT is the **weight** in % of the *total* water in a product. Test methods are eg. the Karl Fischer Titration, or drying balances.


The WATER ACTIVITY,  $a_w$ , is the free water, the only value, which informs about the **micro-biological safety** of a product, specially important for foodstuffs.

In materials it is called the RELATIVE EQUILIBRIUM HUMIDITY (% reh).

The range of 0,00...1,00 $a_w$  corresponds to 0...100% reh.

The relation between moisture content and water activity is the SI, Sorption Isotherm curve, which is *product and temperature specific*. An absolutely constant temperature during all measurements is therefore very important!

Examples of SI-curves:

### **Recording of an SI, Sorption Isotherm curve**

best with a Novasina  $a_w$  instrument with well sealing measuring chamber and full temperature control, selectable from 0°C to 50°C  $\pm 0.2$ K.

Adsorption curve: in order to start, the substance must be dry, so it first has to be totally dried out in an oven (depending on the material eg. 2 hours at 102°C). Then after weighing the dry substance (net weight) on a precision balance, the first possible salt to use is the SAL-T6. After about one hour or more, note the  $a_w$  value and weigh the substance again. When waiting longer, the  $a_w$  value as well as the weight are higher, moving always along the SI-curve, so the “equilibration time” is not critical. Then continue with SAL-T11, up to SAL-T98.

Desorption curve: the substance must be more humid when starting to weigh (to be safely on the desorption curve to avoid SI-hysteresis mistakes), so the first salt to use is SAL-T98. When reaching a value well above 0.9, then SAL-T90 can be placed. After about one hour or more, note the  $a_w$  value and weigh the substance. When waiting longer, the  $a_w$  value as well as the weight are lower, moving always along the SI-curve, so the waiting time is not critical.


After that, place SAL-T75, down to SAL-T6. Then, after totally drying the substance in the oven, the net (dry) weight is known and all points can be plotted.

**How to proceed:**

1. Weigh the tare (empty substance basket and weighing glass)

2. Place the first SAL-T standard

3. Add the SI distance ring

4. Fill basket with substance

5. Place the basket in the ring

6. Close the measuring head

7. Select desired temperature

8. Wait until temperature is reached, then wait for equilibrium (stability factor as desired)

Pictures show "SI set BS" and AW SPRINT. For AWCenters use "SI set aw-Box" and proceed the same way.

**Example** of doing a Sorption Isotherm curve (adsorption) with a **biscuit**:

- Selected temperature at the Novasina instrument: 25,0°C.
- The totally dried out sample has a net weight of 0.485g.
- SAL-T 11% was placed together with the sample in the chamber. After about one hour, the aw-value was 0.087. Weighing the sample showed 0.519g, which is 6.94% more weight; this is the moisture content at this point!
- Next salt: SAL-T 33%. Displayed values: 0.290aw, 0.530g (+ 9.1%).
- Next salt: SAL-T 53%. Displayed values: 0.523aw, 0.553g (+13.84%).
- Next salt: SAL-T 75%. Displayed values: 0.749aw, 0.605g (+24.55%).

Higher humidity was in this case of a "dry biscuit" not of interest.

With these 5 points, an SI-curve can roughly be drawn.

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